

FRONTIER MEDICAL

DURABLE MEDICAL EQUIPMENT AND SUPPLIES

Frontier Main: (801) 979-4638

Fax: (801) 807-0144

Email: orders@ftmed.net

Website: ftmed.net



Table of Contents

Nutrition	pg 3
Respiratory	pg 4-5
Beds and Equipment	pg 6-7
Aid and Mobility	pg 8-9
Incontinence	pg 10
Bathroom Supplies	pg 11
Misc Medical Supplies	pg 12-13
Tips and Maintenance	pg 14-15

NOTE: (*) means the item requires a Letter of Medical Necessity (LMN), Rx, and doctor notes. All other items only need an Rx.

NUTRITION

*Medical Grade Electrolytes



*Nutrition Supplement Drinks



*Food Thickener (powder-based)



*Beneprotein Powder



*Enteral Feeding Pump



*Enteral Feeding Bags



*Gastro Tube (low-pro/std)



Product images for illustration purposes only. Actual product may vary.

RESPIRATORY

* Portable Oxygen Concentrator



* Portable Oxygen Concentrator Backpack



* Stationary Oxygen Concentrator



* Portable Oxygen Concentrator Battery



* Oxygen Tubing



* Nasal Cannulas and Masks



* Pulse Oximeter



* Bubble Humidifier Bottles



* In-Line Water Trap



NOTE: (*) means the item requires an LMN, Rx, and doctor notes. All other items only need an Rx.

RESPIRATORY

* CPAP/BIPAP Machines with Tubing and Masks



Portable Backup Battery (CPAP) with 110 V Outlet



* Nebulizer Machine with Tubing and Masks



* Portable Nebulizer



* Air Purifier



* Room Air Humidifier



Incentive Spirometer



* Peak Flow Device



* Acapella PEP Device



BEDS AND EQUIPMENT

* Hospital Bed Mattress



* Alternating Air Pressure Relief Mattress



* Hospital Bed, Semi or Full-Electric



* Over Bedside Table



Overlays (non-power, dry, gel)



Hospital Bed Rails (half or full length)



* Trapeze Bar



* Tilt-Top Overbed Table



NOTE: (*) means the item requires an LMN, Rx, and doctor notes. All other items only need an Rx.

BEDS AND EQUIPMENT

* Bed Alarm



* Transfer Bed Step



* Bed Foot Support



* Assist-A-Tray



* Large Bed Assist Rail with Organizer



* Small Bed Assist Rail with Organizer



* Bed Cane with Organizer



* Couch/Chair Support Handles



Product images for illustration purposes only. Actual product may vary.

AIDING EQUIPMENT AND MOBILITY DEVICES

Rollator Walker (standard or upright)



*Transport Chair (standard weight or lightweight)



*Mobility Scooter



*Lift Chair



*Wheelchair (standard weight, lightweight, or ultra-lightweight)



Anti-Tip Devices (for wheelchairs)



*Power Wheelchair



*Mobility Ramps



NOTE: (*) means the item requires an LMN, Rx, and doctor notes. All other items only need an Rx.

AIDING EQUIPMENT AND MOBILITY DEVICES

Reacher/Grabber Tool



Cane (single or quad prong)



*Bed Transfer/Positioning Sheet



*Transfer Board



Gait Belt



*Patient Lift (multiple sling options)



*Exercise Bike



*Exercise Treadmill



INCONTINENCE

Briefs, Pull-Up or Tab



Personal Care Wipes



Disposable Exam Gloves



Sani-Cloth Germicidal Wipes



Absorbent Underwear Guards for Men and Women



Disposable Absorbent Bed Pads (Chux)



Catheter Supplies



Male Urinal



Ostomy Supplies



NOTE: (*) means the item requires an LMN, Rx, and doctor notes. All other items only need an Rx.

BATHROOM SUPPLIES

Adjustable Shower Stool (handles or no handles)



Shower Chair (handles or no handles)



Transfer Bench



Sliding Transfer Bench



Commode Chair



Toilet Seat Riser (handles or no handles)



Safety Rails for Toilet



* Grab Bars



* Hand-Held Showerhead



MISC MEDICAL SUPPLIES

*Thermometer
(temporal or oral)



Blood Pressure Monitor



*Weight Scale



Wedge Pillows



Heating Pad



*Tens Unit and accessories



*Knee or Back Brace



Compression (or Diabetic)
Socks and Sleeves



Paraffin Bath (hand/foot wax unit)



*AeroChamber Spacer for Inhaler



*Pneumatic Leg
Compression



*AED Device



NOTE: (*) means the item requires
an LMN, Rx, and doctor notes.
All other items only need an Rx.

MISC MEDICAL SUPPLIES

General Wound Care

Dressings



Split Drain Sponge



Wound Wash



Medical Tape



Hydrocolloid Dressing 4x4



*Tracheostomy Care



*Abdominal Binder



*Truss or Hernia Belt



Tips to Stay Safe at Home

Assess mobility needs in and around the home

- Remove trip hazards.
 - a. De-clutter, remove loose carpeting/rugs or secure with double sided rug tape.
 - b. Mark any steps or changes in floor level with tape or paint in a high-contrast color.
 - c. Remove or secure any electrical or cable cords that might be tripped on.
- Are stair rails secure, weight bearing?
- Assess the need for ramps and stair-lifts.
- Do any doorways need to be widened for wheelchair/walker?

Bathroom safety and accessibility

- Does the shower/tub area need grab bars installed?
- Would they benefit from a shower chair or transfer bench?
- Does the toilet need a surround rail or riser?
- Turn water heater down to reduce scald risk.

Communication

- Would a doorbell camera or intercom be useful?
- Do they have a way to keep a phone with them at all times?
- Would a personal alert system be appropriate?

General safety and convenience considerations

- Store items within reach, do they need a reacher/grabbing device?
- Smoke and carbon monoxide detectors, set a maintenance schedule to check and change batteries.
- Ensure home has adequate lighting.
- Do any areas need a non-slip mat?
- Replace standard doorknobs with lever handles.
- Do any stairways need non-skid treads?
- Ask your Home Health Care company if they can perform a Home Safety Assessment.

The information in this booklet is general in nature; it is not intended as medical advice or a substitute for obtaining medical advice. Frontier Medical makes no representation for the accuracy of information contained herein. The information provided does not override the manufacturers recommendations or directions to maintain the equipment as described in the manufacturer's directions for use.

Basic Oxygen Concentrator and Supply Maintenance

Clean the outside of the oxygen concentrator.

- It is recommended to clean the outside of the concentrator weekly.
- Begin by unplugging the oxygen concentrator from its power source.
- Vacuum vents to prevent air flow obstruction.
- Wipe down the concentrator with a soft damp cloth in a solution of mild dish washing soap and warm water.
- Let the concentrator air-dry or dry with a lint free cloth before plugging into a power source.

Stationary and Portable Oxygen Concentrator Filters

- External Particle Filters: Clean your external particle filters weekly by rinsing them with warm water and let them air dry.
- Internal Filters: You can request internal filters for the Stationary and Portable every two years.

Suggested Timeline for Changing O2 Supplies

Nasal cannulas/masks

- Replace weekly.

Ear cushions and swivel connectors

- Replace monthly.

Water traps

- Every day, or as needed, twist/pull water trap tube open to empty it.
- The placement of the water trap should be next to the nasal cannula not near the humidifier bottle.
- Replace monthly.

Humidifier bottle

- Only use distilled water.
- Weekly wash the humidifier bottle with soap and warm water.
- Replace humidifier bottle monthly.

Green extension tubing

- Replace every other month.

Additional Tips

- Do NOT submerge the concentrator in water.
- Do not place liquids near your concentrator.
- Keep concentrator 6-12 inches away from walls, furniture, and curtains to maintain airflow.
- Keep away from heat or open flame.
- Don't smoke or allow others to smoke near an oxygen concentrator.

Because Your Needs Matter.

We work hard to help you stay
at home, through quality products
and compassionate customer service.
Let us know how we can help!

The offerings in this booklet and images are just examples of common
equipment we can provide. We are not limited to what is listed here and
we can provide nearly all home medical equipment to meet your specific needs.

Frontier Main: (801) 979-4638

Fax: (801) 807-0144

Email: orders@ftmed.net

Website: ftmed.net