

Frontier Medical

Durable Medical Equipment & Supplies



*"Because your
needs matter!"*

Main: (801) 979-4638

Fax: (801) 807-0144

Email: orders@ftmed.net

Website: www.ftmed.net

Visit our new website!

See our patient resources,
online catalog & more.



Our Mission


At Frontier Medical, we deeply value our relationships and appreciate your choice to trust us among other options for durable home medical equipment. We are committed to personalized service and high-quality supplies to enhance your home life.

We specialize in assisting DEEOIC, DFEC, and the DOL's OWCP Black Lung Program patients, with the ability to work directly with DOL-approved physicians and preferred Home Health Care providers across the United States.

Products Designed for Your Health

We provide a diverse selection of reliable, user-friendly products from well-known national brands. We are dedicated to ensuring you receive the essentials for your health - delivered straight to your door at no cost to you.

Where We Serve



** At this time we cannot serve patients in Alabama and Maryland.
Subject to change.

We provide top-of-the-line home medical supplies to those who qualify. Our knowledgeable staff is eager to help you get the home medical equipment you need as soon as possible.

With thousands of clients nationwide and projected growth, we are here to help in any way we can - because your needs matter.



Our Clients Say...

"Anytime I ask you guys for anything I usually get my things in a few days, it's such a fast turnaround."

~ Client, December 2025

"Thank you again for everything. God has truly blessed my wife and I through all of you at Frontier."

~ Client, November 2025

"Thank you all for being the most wonderful DME company I've ever had the pleasure to work with!"

~ Client, July 2023

"I cannot say enough nice things about the Frontier Medical team! Everyone there is so professional, kind, patient, etc."

~ Client, June 2022

"I have never worked with so many incredible people dedicated to detail and helping others. Good job everyone!"

~ Client, April 2024



** To ensure HIPAA compliance, our reviews are censored to protect the identities of our clients while providing you with honest feedback.

Table of Contents

Nutrition	pg 5
Respiratory	pg 6-9
Beds, Equipment & Safety	pg 10-11
Aid & Mobility	pg 12-14
Bathroom Supplies	pg 15
Personal Care Items	pg 16
Misc Medical Supplies	pg 17-18
DME Checklist	pg 19

NUTRITION



Medical Grade Electrolytes



Medical Grade Electrolytes



Medical Grade Electrolytes



Nutrition Supplement Drinks



Nutrition Supplement Drinks



Nutrition Supplement Drinks



Food Thickener (Powder-Based)



Beneprotein Powder



Enteral Feeding Pump



Enteral Feeding Bags



Gastro Tube (Low-Pro/Standard)

RESPIRATORY



*

Stationary Oxygen Concentrator



*

Portable Oxygen Concentrator



*

Portable Concentrator Accessories
(Backpack, Battery, Power Cords, etc.)



*

Oxygen Supplies
(Nasal Cannulas, Oxygen Tubing, etc.)



*

Acapella PEP Device



*

Peak Flow Device



Incentive Spirometer




*

Pulse Oximeter



*

Room Humidifier



*

Air Purifier



**Don't See What
You're Looking For?**

Please feel free to call and speak to one of our representatives! We are not limited to what is listed here and we can provide nearly all home medical equipment to meet your specific needs.

**Call Us Today at
(801) 979-4638**

Basic Oxygen Concentrator and Supply Maintenance

Cleaning the Oxygen Concentrator:

- Clean the exterior of the oxygen concentrator once per week.
- Before cleaning, unplug the concentrator from the power source.
- Use a vacuum to remove dust from the vents to prevent airflow obstruction.
- Wipe the exterior with a soft, damp cloth using a solution of mild dishwashing soap and warm water.
- Allow the concentrator to air dry, or dry it with a lint-free cloth, before plugging it back in.

Stationary and Portable Oxygen Concentrator Filters:

- External Particle Filters - Clean weekly by rinsing with warm water. Allow them to air dry completely before reinstalling.
- Internal Filters - Internal filters for both stationary and portable units should be changed every two years and can be requested as needed.

Suggested Timeline for Changing Oxygen Supplies

Nasal Cannulas / Masks

- Replace weekly.

Ear Cushions, Swivel Connectors, and Humidifier Tubing

- Replace monthly.

Water Traps

- Empty daily, or as needed, by twisting or pulling open the water trap tube.
- Position the water trap near the nasal cannula, not near the humidifier bottle.
- Replace monthly.

Humidifier Bottle

- Only use distilled water.
- Ensure the lid is properly threaded to prevent leaks.
- Wash weekly with soap and warm water.
- Replace monthly.

Green Oxygen Tubing

- Replace every other month.

Additional Tips

- Do not submerge the concentrator in water.
- Keep liquids away from the concentrator.
- Maintain at least 6-12 inches of clearance from walls, furniture, and curtains to allow proper airflow.
- Keep away from heat sources or open flames.
- Do not smoke or allow others to smoke near the oxygen concentrator.

Basic CPAP / BiPAP and Supply Maintenance

Caring for Your CPAP / BiPAP Supplies

To ensure proper use of your PAP Therapy, it's important to regularly inspect, clean and care for your equipment.

Mask Care

- Facial oils can accumulate on your mask and reduce the performance of the seal. It is recommended to wash your face with water to remove any excess facial oils before wearing your mask each night.

Device Care

- Wipe the device down and clean its air tubing weekly. It is recommended to check the air filter at least every few weeks for dust build-up. If the air filter is blocked by dust, it can reduce the air flow and diminish air pressure being delivered to you.

Humidifier Care

- Wash the humidifier tub in warm water using mild detergent (like dish soap) and allow to air dry. Inspect the water tub once a month for wear and deterioration as a worn out tub can cause leaking. Always replace the tub if it is leaking, cracked or has become cloudy.

** Please refer to the relevant product user guide provided with your device and supplies for more information.

Important

Only ever wash supplies by hand with a mild detergent to avoid damage and harmful residues.

A list of common cleaning solutions to avoid:

- × Aromatic-based solutions or scented oils
- × Antibacterial soaps
- × Alcohol
- × Bleach
- × Strong-smelling products

Replacing Your CPAP / BiPAP Supplies

Every 2 Weeks

- Filters
 - Replace more often if there are any holes or blockages by dirt or dust.
- Nasal Pillows / Cushions

Monthly

- Full Face Mask Cushion

Every 3 Months

- Air Tubing / Hose
- Mask Frame
- Mask Headgear
- Humidifier Tub

Talk to one of our CPAP representatives to sign up for a recurring order of supplies!
Call Us Today at (801) 979-4638

RESPIRATORY



CPAP / BiPAP Machine
with Tubing & Masks



AirMini Travel CPAP
(Travel Case, Hose System, HumidX, etc.)



CPAP / BiPAP Supplies
(Mask, Hose, Water Chamber, Filters, etc.)



Portable Backup Battery
110 V Outlet (TSA & FAA Approved)



Portable Backup Battery
120 V Outlet (12-16 hr charge)



CPAP Therapy Pillow



Nebulizer Machine



Portable Nebulizer Machine



Nebulizer Supplies
(Mouthpiece, Tubing & Mask)

BEDS, EQUIPMENT & SAFETY



Hospital Bed
(Semi or Full-Electric)



Hospital Bed Mattress



Alternating Air Pressure
Relief Mattress



Mattress Overlay
(Non-Power, Dry or Gel)



Hospital Bed Rails
(Half or Full Length)



Transfer Bed Sheet



Bed Cradle & Foot Support



Transfer Bed Step



Couch / Chair Support Handles



Over Bedside Table



Tilt-Top Over Bedside Table



Assist-A-Tray



Bed Alarm



Pull Cord Alarm



Home Emergency Medical
Alert System (Wearable)

Tips to Stay Safe at Home

Assess Mobility Needs in and Around the Home:

- Remove trip hazards
 - a. De-clutter, remove loose carpeting / rugs or secure with double sided rug tape.
 - b. Mark any steps or changes in floor level with tape or paint in a high-contrast color.
 - c. Remove or secure any electrical or cable cords.
- Are stair rails secure and/or weight bearing?
- Assess the need for ramps and stair-lifts.
- Do any doorways need to be widened for wheelchairs / walkers?

Bathroom Safety and Accessibility:

- Does the shower / tub area need grab bars installed?
- Would they benefit from a shower chair or transfer bench?
- Does the toilet need a surround rail or riser?
- Turn water heater down to reduce scald risk.

Communication:

- Would a doorbell camera or intercom be useful?
- Do they have a way to keep a phone with them at all times?
- Would a personal alert system be appropriate?

General Safety and Convenience Considerations:

- Store items within reach, do they need a reacher / grabbing device?
- Smoke and carbon monoxide detectors, set a maintenance schedule to check and change batteries.
- Ensure home has adequate lighting.
- Do any areas need a non-slip mat?
- Replace standard doorknobs with lever handles.
- Do any stairways need non-skid treads?
- Ask your Home Health Care company if they can perform a Home Safety Assessment.

The information in this booklet is general in nature; it is not intended as medical advice or a substitute for obtaining medical advice. Frontier Medical makes no representation for the accuracy of information contained herein. The information provided does not override the manufacturers recommendations or directions to maintain the equipment as described in the manufacturer's directions for use.

AID & MOBILITY



*

Trapeze Bar



*

Patient Lift
(Multiple Sling Options)



*

Transfer Board



Gait Belt



Freestanding Cane



Quad-Based Cane



Crutches



Reacher / Grabber Tool



*

Mobility Ramps



Rollator Walker



Folding Walker



Upright Walker

AID & MOBILITY



*

Wheelchair



*

Portable Power Wheelchair



*

Power Wheelchair



*

Wheelchair Seat / Back Pads



*

Wheelchair Arm Rest Pads



*

Wheelchair Elevating Leg Rests



*

Wheelchair Anti-Tippers



*

Wheelchair Safety Belts



*

Wheel Lock Extensions



*

Transport Chair



*

Lightweight 3-Wheel Scooter



*

Mobility Scooter


AID & MOBILITY



Pedal Exerciser




Rowing Machine



Elliptical Machine



Exercise Treadmill



Exercise Bike



Lift Chair



Don't See What You're Looking For?

Please feel free to call and speak to one of our representatives! We are not limited to what is listed here and we can provide nearly all home medical equipment to meet your specific needs.

**Call Us Today at
(801) 979-4638**

BATHROOM SUPPLIES



Adjustable Shower Stool
(Handles or No Handles)



Shower Chair
(Handles or No Handles)



Transfer Bench



Sliding Transfer Bench



Hand-Held Showerhead



Grab Bars



Safety Rails for Toilet



Toilet Seat Riser
(Handles or No Handles)



Commode Chair

PERSONAL CARE ITEMS



Ostomy Supplies



Catheter Supplies



Disposable Absorbent Bed Pads (Chux)



Absorbent Underwear Guards



Briefs (Pull-Up or Tab)



Panty Liners



Disposable Exam Gloves



Personal Care Wipes



Sani-Cloth Germicidal Wipes



General Wound Care

MISC. MEDICAL SUPPLIES



Back Brace



Knee Brace



Compression / Diabetic Sleeves



Compression / Diabetic Socks



Pneumatic Leg Compression



Hot / Cold Water Therapy



Truss or Hernia Belt



Abdominal Binder



AED Device



Tracheostomy Care



Suction Machine & Supplies



TENS Unit & Accessories


MISC. MEDICAL SUPPLIES



Inhaler Spacer



Heating Pad



Weight Scale



Blood Pressure Monitor



Temporal Thermometer



Oral Thermometer



Anti-Bedsore Pillow



Wedge Pillow



Knee Pillow



Half-Moon Pillow

DME Checklist

Personal Care / Wound Care

- ☐ Incontinence Garments
- ☐ Chux / Absorbent Pads
- ☐ Gloves
- ☐ Perineal Spray
- ☐ Zinc Oxide
- ☐ Personal Care Wipes
- ☐ Germicidal Wipes
- ☐ Alcohol Wipes
- ☐ Urinal
- ☐ Catheter
- ☐ Bed Pan
- ☐ Wound Care
- ☐ Ostomy Supplies

Bathroom Safety

- ☐ Toilet Seat Riser
- ☐ Toilet Safety Rail
- ☐ Commode Chair
- ☐ Non-Slip Bathmat
- ☐ Shower Stool
- ☐ Shower Chair
- ☐ Shower Bench
- ☐ Sliding Shower Bench
- ☐ Grab Bars
- ☐ Hand-Held Shower Head

Nutrition

- ☐ Ensure
- ☐ Boost
- ☐ Electrolytes
- ☐ Glucerna
- ☐ Food Thickener
- ☐ Beneprotein Powder
- ☐ Enteral Feeding Pump

Vitals Equipment

- ☐ Blood Pressure Monitor
- ☐ Stethoscope
- ☐ Pulse Oximeter
- ☐ Temporal Thermometer
- ☐ Oral Thermometer
- ☐ Scale

Oxygen

- ☐ Stationary O2 Machine
- ☐ Portable O2 Machine
- ☐ Portable O2 Charger
- ☐ Portable O2 Battery
- ☐ Portable O2 Backpack
- ☐ Nasal Cannulas
- ☐ Oxygen Tubing
- ☐ Humidifier Bottle
- ☐ Misc. O2 Supplies
- ☐ Stationary O2 Filters
- ☐ Portable O2 Filters

Beds

- ☐ Hospital Bed (Full / Semi Electric)
- ☐ Hospital Bed Mattress
- ☐ Hospital Bed Rails (Full / Half)
- ☐ Trapeze Bar
- ☐ Overbed Table (Standard / Tilt)
- ☐ Alternating Pressure Mattress
- ☐ Alternating Pressure Pad
- ☐ Bed Alarm
- ☐ Adjustable Bed
- ☐ Gel Pressure Mattress
- ☐ Bed Sheets
- ☐ Mattress Protector
- ☐ Bed Support
- ☐ Bed Stepstool
- ☐ Couch Support Handles

Mobility Devices

- ☐ Walker (Standard / Rollator / Upright)
- ☐ Cane / Quad Cane
- ☐ Crutches
- ☐ Lift Chair
- ☐ Patient Lift
- ☐ Wheelchair
- ☐ Power Wheelchair
- ☐ Scooter
- ☐ Transport Chair
- ☐ Scooter Lift
- ☐ Ramp

Misc. Items

- ☐ Wedge / Positioning Pillow
- ☐ Back / Seat Cushion
- ☐ Heating Pad
- ☐ Reacher/Grabber
- ☐ TENS Unit & Supplies
- ☐ Sheep Skin Pad
- ☐ Compression Socks / Sleeves
- ☐ Non-Slip Socks
- ☐ Diabetic Socks
- ☐ Transfer Board / Sheet
- ☐ Gait Belt
- ☐ AED Device
- ☐ Suction Machine
- ☐ Assist-a-Tray
- ☐ Pneumatic Compression
- ☐ Hernia Belt
- ☐ Braces
- ☐ Abdominal Binder
- ☐ Exercise Equipment
- ☐ Tracheostomy Supplies

Misc. Respiratory Items

- ☐ Spirometer
- ☐ Acapella Device
- ☐ Peak Flow Device
- ☐ Air Purifier & Filters
- ☐ Humidifier
- ☐ Inhaler Spacer

Nebulizer

- ☐ Desktop Nebulizer
- ☐ Portable Nebulizer
- ☐ Nebulizer Supplies

CPAP

- ☐ CPAP Machine
- ☐ BiPAP Machine
- ☐ Travel CPAP Machine
- ☐ CPAP / BiPAP Supplies
- ☐ Portable Backup Battery
- ☐ CPAP Pillow

COMING SOON

Watch for new kit options on our Website!

Because Your Needs Matter.

We work hard to help you stay safe
at home, through quality products
and compassionate customer service.

Let us know how we can help! (801) 979-4638



The offerings in this booklet and images are just examples of common equipment we can provide. We are not limited to what is listed here and we can provide nearly all home medical equipment to meet your specific needs.

Main: (801) 979-4638
Fax: (801) 807-0144
Email: orders@ftmed.net
Website: www.ftmed.net

Visit our new website!

See our patient resources,
online catalog & more.

